

# Analysis of Halogenated Acetic Acids and Dalapon in Drinking Water Using EPA Method 552.1 with Semi-Automated Solid Phase Extraction (EZSpe®)

## Introduction

Haloacetic acids are known by-products of water chlorination. They often accumulate in significant quantities and have been implicated in multiple undesirable health outcomes by various epidemiological analyses.

EPA Method 552.1 outlines the procedure for the extraction and analysis of these compounds in drinking water. The extraction method outlines the use of solid phase extraction for water matrix samples employing both cartridges and disks. Consistent with other EPA 500 series methods, EPA 552.1 incorporates a rigid set of QC and acceptance criteria requiring precise and reproducible analytical practices. The potential for error and the variability associated with manual extractions makes the benefits of semi-automating these processes apparent.

To meet demands for a low cost method that requires less financial investment than fully automated systems, FMS developed a simple semi - automated system which is fast, inexpensive and yields high quality data.

## Instrumentation

- FMS EZSpe® System
- FMS SuperVap®
- Vacuum pump
- Agilent 7890A GC FID/ECD

## Consumables

- FMS, Inc. 1 g SDVB cartridge
- FMS sodium sulfate cartridge
- Ultra pure DI water
- Fisher 6 N Hydrochloric Acid/Sodium Hydroxide
- Fisher 6 N sulfuric acid
- Fisher Pesticide Grade Methanol
- Fisher methyl tert-butyl ether
- Restek 552.1 spiking standards

## Procedure

- 6 samples (100 mL water each) are prepared and acidified with HCl till pH ~ 5
- Spike with various 552.1 standards
- Put sample bottles in place and fill rinse bottles with 4 mL 10% sulfuric acid in methanol
- Cartridges are installed in each of the six positions.

### Stage 1:

- Vacuum is turned on
- Cartridges are conditioned with 10 mL methanol, 10 mL water, 10 mL 1M HCl in methanol, 10 mL water, 10 mL 1M NaOH, 10 mL water (keep wet after each addition)
- Samples are loaded across cartridges under vacuum, at 2 mL/min.
- Cartridges are dried under vacuum with 10 mL methanol
- Sample bottles are automatically rinsed from the rinse bottles with 4 mL 10% sulfuric acid in methanol

### Stage 2:

- Sulfuric acid/methanol from sample bottles is loaded across the cartridges and
- To the eluent is added 2.5 mL of methyl tert-butyl ether, which is then heated for 1 hour and collected for analysis into Direct to GC Vial Collection Vessels
- Extracts are dried over sodium sulfate or in line cartridges can be used downstream from SDVB cartridges

### FMS SuperVap®

- Pre-heat temp: 45 °C
- Pre-heat time: 15 minutes
- Heat in Sensor mode at 45 °C under nitrogen (7-10 psi)
- Direct to GC Vial Vessel Reduce to 1 mL
- Samples are now ready for analysis




Table 1 with recoveries for the 552.1 compounds

Compound	%Recovery	Stdev
Monochloroacetic acid	83.2	1.63
Dichloroacetic acid	91.4	2.75
Trichloroacetic acid	76.8	1.99
Monobromoacetic acid	88.2	1.92
Bromochloroacetic acid	96.6	1.84
Dibromoacetic acid	100.3	2.03
Dalapon	79.2	3.60

### Conclusions

Reviewing the sample data shows high recoveries for the seven 552.1 analytes, demonstrating excellent efficiency for these compounds. Samples can be taken from collection bottle to GC vial in one quick, consistent, reproducible process that will save laboratories both time and money.



For more information contact FMS:  
FMS, Inc.  
580 Pleasant Street  
Watertown, MA 02472  
Phone: (617) 393-2396  
Fax: (617) 393-0194  
Email: [onlineinfo@fms-inc.com](mailto:onlineinfo@fms-inc.com)  
Web site: [www.fms-inc.com](http://www.fms-inc.com)

FMS EZSpe® System