


Optimized Sample to Result Workflow of Dioxin and PCBs Automated Sample Prep and GC-MS/MS Analysis in Environmental and Food

Reinvigorating Dioxins Analysis with GC/MS/MS


Joel Ferrer
Product Manager, Triple Quadrupole GC/MS




Introducing the 7000E GC/TQ and 7010C GC/TQ



7000E GC/TQ



7010C GC/TQ




Next-Generation Mass Spec Intelligence

Intelligence Powered Advances on GC/TQ



NEW SWARM Autotune completes 2X faster





System Screening & System Evaluation assess presence of ions and electronics status



Diagnostic Tune used for troubleshooting



A detailed system report (MSR) can be sent to a service engineer for speedy diagnostic efforts



Introduction

Persistent Organic Pollutants (POPs) in Stockholm Convention


- The Stockholm Convention on POPs (2001) is a global treaty to protect **human health** from chemicals that remain in the environment and are **persistent, bioaccumulative** and transportable across the globe.
- Include:
 - Industrial chemicals ex. PCBs, hexachlorobenzene
 - Pesticides ex. Aldrin, DDT, endrin, toxaphene etc.
 - Pharmaceuticals
 - Solvents
 - By-products ex. **Dioxins & Furans**
- Initially 'dirty dozen' but new POPs include Perfluoro alkyl substances (**PFAS**), short chain chlorinated paraffins (**SCCPs**), pesticides, Deca-BDE etc.
- Provisions include eliminating production and intentional uses, managing and disposing wastes in a suitable environmentally safe manner



A History of Dioxin Analysis and USEPA 1613B

Most Widely Used Dioxin Method Globally

- A **HRGC/HRMS** method that is used to analyze 2,3,7,8-tetrachloro dibenzo-p-dioxin (TCDD) in municipal and industrial discharges
- Used to determine **PCDDs and PCDFs** in water, soils, sludges and other matrices
- Specificity is provided for determination of the **seventeen 2,3,7,8-substituted PCDDs and PCDFs**
- This method is also used for the analysis of **drinking water, food products, and human/animal tissue** samples
- **Many regions** outside the U.S. also use the method as guidance or strictly follow it




PCDD and PCDF Nomenclature and Isomers

EPA 1613B Analytes

Dioxins

Polychlorinated dibenzodioxin (PCDD)




2,3,7,8 – TCDD

(2,3,7,8-Tetrachlorinated dibenzodioxin)

Furans

Polychlorinated dibenzofurans (PCDF)



2,3,4,7,8 – PCDF

(2,3,4,7,8-Pentachlorinated dibenzofurans)

Chlorine Atoms	PCDD Isomers	PCDF Isomers
1	2	4
2	10	16
3	14	28
4	22	38
5	14	28
6	10	16
7	2	4
8	1	1
Total	75	135

75 total PCDD isomers

135 total PCDF isomers

17 toxic isomer

Toxic Equivalent Factors (TEF)

17 Toxic Dioxin Compounds

Toxic PCDDs	I-TEF	WHO ₂₀₀₅ -TEF
2,3,7,8-TetraCDD	1	1
1,2,3,7,8-PentaCDD	0.5	1
1,2,3,4,7,8-HexaCDD	0.1	0.1
1,2,3,6,7,8-HexaCDD	0.1	0.1
1,2,3,7,8,9-HexaCDD	0.1	0.1
1,2,3,4,6,7,8-HeptaCDD	0.01	0.01
1,2,3,4,6,7,8,9-OctaCDD	0.001	0.0003


Toxic PCDFs	I-TEF	WHO ₂₀₀₅ -TEF
2,3,7,8-TetraCDF	0.1	0.1
1,2,3,7,8-PentaCDF	0.05	0.03
2,3,4,7,8-PentaCDF	0.5	0.3
1,2,3,4,7,8-HexaCDF	0.1	0.1
1,2,3,6,7,8-HexaCDF	0.1	0.1
1,2,3,7,8,9-HexaCDF	0.1	0.1
2,3,4,6,7,8-HexaCDF	0.1	0.1
1,2,3,4,6,7,8-HeptaCDF	0.01	0.01
1,2,3,4,7,8,9-HeptaCDF	0.01	0.01
1,2,3,4,6,7,8,9-OctaCDF	0.001	0.0003

$$TEQ = \sum_{n=1}^{17} ([PCDD/F]_i (ng/L) \times TEF_i) (ng TEQ/L)$$


Environmental Dioxins Testing Today

Agilent 7010C Triple Quadrupole GC/MS

- The U.S. EPA has approved an alternate testing protocol for 1613B – Analysis of Dioxins
- The new alternate testing protocol accepts GC/TQ as an **equivalent technology** to analyze dioxins in Environmental Samples



**U.S. EPA
1613B**



**SGS AXYS
16130**

Comparison of U.S. EPA 1613B vs SGS AXYS 16130

Requirement		U.S. EPA 1613B	SGS AXYS 16130
No Change	Valid for EPA compliance reporting	Yes	Yes
	Analyte List	Tetra through Octa Dioxins and Furans (210 congeners – 17 toxics)	Same as 1613B
	Matrices validated	Non-potable water, soil, sediment, fish tissue and other solids	Same as 1613B
	Sample Prep.	Soxhlet extraction/SPE + optional carbon clean-up	Same as 1613B
	Reporting units	Toxicity equivalents (TEQs)	Same as 1613B
Updated	Method Detection Levels (MDLs)	4.4 pg/L for 2,3,7,8-TCDD	1.1 pg/L for 2,3,7,8-TCDD (all MDLs calculated and in method)
	Detector	GC/HRMS (R>10,000 at 10% valley ex. Magnetic Sector)	GC/TQ (Unit resolution)
	Instrument tuning requirements	PFK tune to meet resolution and m/z deviation of <5 ppm	Standard Autotune to meet manufacturer's specifications of mass accuracy
	Lock mass reference requirement	Infusion of PFTBA	No lock mass; calibrant Infusion of PFTBA 'at manufacturer's specified limits' for matrix evaluation

Dioxins via SGS AXYS 16130

With this approval, some modification is required...


6.10 Mass spectrometer—Tandem quadrupole mass spectrometer equipped with either an atmospheric pressure (AP) ionization source or an electron ionization source (70 eV range) capable of repetitively selectively monitoring at least 32 transitions at unit resolution during a period of approximately 1 second, and shall meet all of the performance specifications in Section 10.

The MS/MS must have a mechanism for constantly bleeding PFTBA into the source during the analytical run, such as a needle valve or a capillary bleed directly into the source enclosure.

10.2.1.2 The response of the MRM transition shown in Table 8 for the reference compound (PFTBA) is monitored throughout the run. The response of the MRM product ion shall not vary by more than $\pm 20\%$ throughout the run. Variations by more than 20% within the possible elution window of a CDD/CDF indicate the presence of coeluting interferences that may significantly alter the response of the mass spectrometer. Reinjection of another aliquot of the sample extract will not resolve the problem. Additional cleanup of the extract may be required to remove the interferences.

System Configuration

Reference Compound Introduction Valve (RCI Valve)




NEW
Reference Compound
Introduction Valve
(RCIV)

**Calibration
Valve**


System Configuration

Also compatible with CTC PAL for sample introduction (1 mL vial)



EPA Method 1613B Sample Preparation

No changes to sample preparation or collection for SGS AXYS 16130



EPA 1613B is a performance-based method

Chromatography Performance




Figure 1. MRM chromatograms for TCDFs, labeled TCDF ISTD, TCDDs, and labeled TCDD ISTD.





Figure 2. MRM chromatograms for HxCDDs and the corresponding ISTD.

Chromatography Performance



Method Table														
Time Segment: <All>			Compound: 2,3,7,8-TCDD			Reset Table View								
	Quantifier													
	Name	TS	Transition	Scan	Type	Resolution Calculation Type		Resolution Limit						
	2,3,7,8-TCDD	1	319.9 -> 256.9...	MRM	Target	Valley Height/Peak Heig...		25.0						
Sample							2,3,7,8-TCDD Results							
Name	Data File	Type	Level	Vial	Acq. Date-Time	Acq. Method File	RT	Resp.	MI	Calc. Conc.	Accuracy	S/N	Resolution F.	Resolution R.
DX041D-CAL_01-73	DX9Z0444.D	Cal	CS3	7	8/22/2019 2:31 AM	TQEI_DB5_DX_11	26.351	221151		9.0183	90.2	2339.21	20.4	7.8

Figure 3. (Left)
2,3,7,8-TCDD and its close eluters.

Figure 4. (Above)

Top Row: Method setup for resolution check in MassHunter Quantitative Analysis.

Bottom Row: Front and rear valley height/peak height resolution calculated for 2,3,7,8-TCDD and its closest eluting isomers.

Reference Compound Stability

The need for lock mass monitoring of the GC/HRMS system for Method 1613B was replaced by use of a stability reference compound in the GC/TQ method.

The Agilent Reference Compound Introduction Valve (RCI Valve):

- Provides an **optimized flow** of PFTBA to observe any changes in the ionization efficiency and ion transmission
 - Seen as a change in the PFTBA signal intensity
- **No software changes or updates** are required– Controlled easily through Agilent MassHunter software




Figure 5. Stability of PFTBA response over the entire analytical run.

Compatible on 7010 Series GC/TQ

Method detection Limits (MDLs) with 7890/7010 GC/TQ Far Surpassed Method 1613B MDLs

The GC/TQ MDL results for the aqueous (1 L), solid (10 g), and tissue (10 g) samples:

Compound	Aqueous	Solid	Tissue
	MDL and (MRL) in pg/L	MDL and (MRL) in pg/g	MDL and (MRL) in pg/g
2,3,7,8-TCDD	1.1 (10)*	0.029 (1)	0.057 (0.5)
1,2,3,7,8-PeCDD	1.39 (50)	0.037 (5)	0.051 (2.5)
1,2,3,4,7,8-HxCDD	1.05 (50)	0.042 (5)	0.061 (2.5)
1,2,3,6,7,8-HxCDD	1.08 (50)	0.045 (5)	0.033 (2.5)
1,2,3,7,8,9-HxCDD	1.78 (50)	0.064 (5)	0.067 (2.5)
1,2,3,4,6,7,8-HpCDD	1.19 (50)	0.070 (5)	0.032 (2.5)
OCDD	9.4 (100)	0.311 (10)	0.085 (5)
2,3,7,8-TCDF	0.56 (10)	0.60 (1)	0.056 (0.5)
1,2,3,7,8-PeCDF	1.0 (50)	0.037 (5)	0.046 (2.5)
2,3,4,7,8-PeCDF	1.25 (50)	0.039 (5)	0.033 (2.5)
1,2,3,4,7,8-HxCDF	0.89 (50)	0.032 (5)	0.029 (2.5)
1,2,3,6,7,8-HxCDF	1.11 (50)	0.031 (5)	0.046 (2.5)
1,2,3,7,8,9-HxCDF	1.22 (50)	0.048 (5)	0.084 (2.5)
2,3,4,6,7,8-HxCDF	1.26 (50)	0.026 (5)	0.034 (2.5)
1,2,3,4,6,7,8-HpCDF	0.92 (50)	0.255 (5)	0.064 (2.5)
1,2,3,4,7,8,9-HpCDF	1.35 (50)	0.028 (5)	0.043 (2.5)
OCDF	2.81 (100)	0.365 (10)	0.113 (5)


*The current MDL for 2,3,7,8-TCDD with the GC/HRMS 1613B is 4.4 pg/L

The MDL achieved with Agilent 7010B GC/TQ is **4 times** better!

Total PCDDs and PCDFs Reporting

Total PCDD and PCDF concentrations from the real-world sample extracts were reported by MassHunter software for each level of chlorination by **summing the concentration of the individual peaks meeting quantification criteria** (peak shape, $S/N \geq 2.5$, and product ion ratio $\pm 10\%$) in the appropriate retention time window.

Comparable total PCDDs and PCDFs reported using **GC/HRMS** and **GC/TQ**:



Conclusions

GC/TQ provides many of the specificity and sensitivity advantages of HRMS

Low Cost of Ownership

Lower instrument cost without the need for expensive consumables, maintenance, or specialized instrument operators

Versatility

7010C can be used strategically with other EPA methods leading to faster ROI

Flexibility

7010C utilizes electron ionization meaning it is compatible with several established libraries including the NIST GC Library



For More Information...

Application Note

An Alternate Testing Protocol for EPA 1613B using Agilent Triple Quadrupole GC/MS

Determination of 2,3,7,8-substituted tetra- through octa-chlorinated dibenzo-p-dioxins and dibenzofurans

Author

Coreen Hamilton and Xinhui Xie,
SGS AXYS Analytical Services Ltd.

Tarun Anumol,
Anastasia Andrianova, and
Dale Walker,
Agilent Technologies, Inc.



Abstract

This study provides data used to create an alternate testing protocol for the U.S. Environmental Protection Agency (EPA) to use for Agilent 7010B Triple Quadrupole GC/MS analysis of tetra- through octa-dioxins and furans that is equivalent to EPA Method 1613B. EPA Method 1613B is used for the determination of the 17 toxic tetra- through octa-chlorinated Dibenzo-p-Dioxins and Dibenzofurans (CDDs/CDFs) in aqueous, solid, and tissue matrices by isotope dilution gas chromatography/high-resolution mass spectrometry (GC/HRMS) using magnetic sector instruments. Traditionally used for dioxins analysis because of their high sensitivity, GC/HRMS instruments are expensive to maintain, require a highly specialized skill set to operate, and are being phased out by manufacturers.

Agilent
CrossLab
From Insight to Outcome

Analysis of Dioxins in Environmental Samples using GC/MS

Consumable Workflow Ordering Guide



Visit us at www.Agilent.com