

Automated Solid Phase Extraction of 40 native PFAS Compounds in Wastewater Using Vacuum and Positive Pressure

Tom Hall

thall@fms-inc.com

Fluid Management Systems

Billerica, MA USA

Purpose

- Sample Preparation Workflows and tools for EPA method 1633 (draft)
- Comparison of Automated and Semi-Automated Solid Phase Extraction for PFAS applications
- EPA Method 1633 applications



Target Analytes

Method	533	537.1	1633	ISO21675
PFBA	x		x	x
PFPeA	x		x	x
PFHxA	x	x	x	x
PFHpA	x	x	x	x
PFOA	x	x	x	x
PFNA	x	x	x	x
PFDA	x	x	x	x
PFUnA	x	x	x	x
PFDoA	x	x	x	x
PFTTrDA		x	x	x
PFTeDA		x	x	x
PFBS	x	x	x	x
PFPeS	x		x	
PFHxS	x	x	x	x

Target Analytes

Method	533	537.1	1633	ISO21675
PFHpS	X		X	X
PFOS	X	X	X	X
PFNS			X	
PFDS			X	X
PFDoS			X	
4:2FTS	X		X	
6:2FTS	X		X	X
8:2FTS	X		X	X
PFOSA			X	X
NMeFOSA			X	X
NEtFOSA			X	X
NMeFOSAA		X	X	X
NEtFOSAA		X	X	X
NMeFOSE			X	

Target Analytes

Method	533	537.1	1633	ISO21675
NETFOSE			X	
HFPO-DA	X	X	X	X
ADONA	X	X	X	X
PFMPA	X		X	
PFMBA	X		X	
NFDHA	X		X	
9CI-PF3ONS	X	X	X	X
11CI-PF3OUdS	X	X	X	
PFEESA	X		X	
3:3 FTCA			X	
5:3 FTCA			X	
7:3 FTCA			X	
PFOcDA				X
PFHxDA				X
8:2 FTUCA				X
8:2 diPAP				X

Optimizing the PFAS Analysis Workflow

- **Automate the Sample Prep Workflow**
 - Automate the Solid Phase Extraction Step
 - Automate the Concentration/Evaporation Step
- **Automated, Semi-Automated SPE extractions and Concentration**
 - Reduces Human Error
 - Reduces Outside contamination and background
 - Reduces Solvent Usage
 - Reduces Labor
- **Use SPE solutions to deliver consistent, reproducible results**
- **Different configurations to handle Matrix type and Budget**
- **No Teflon Components**

Comparison of Manual SPE vs. Automated SPE Methods

Manual SPE

Manually Separate Waste

<10 mls solvent evaporate

Run times are ~ 30 to 90 minutes

Technician Time 25 minutes

Physical transfer and Concentration steps

Concentration steps 70 to 200 minutes

Automated/Semi-Auto SPE

Separates Aqueous and Organic Waste

<10 mls solvent to evaporate

Run times are ~20 to 50 minutes or less

Technician time 5 minutes

Automatic and Direct to Concentration
delivery and completion

Concentration step 70 minutes ready for
injection



Determining Factors

- **Closed System**
- **Ability to load samples by vacuum or positive pressure consistently.**
- **Handle difficult Matrices**
- **Ability to dry cartridges by vacuum and positive Nitrogen gas pressure.**
- **Easily handle various cartridge designs and sizes without cumbersome modifications.**
- **Easy, Simple Sample delivery**
- **Bottle Rinse**



Automated Solid Phase Extraction front end for LC/MS



EconoTrace® PFC



TurboTrace® PFC



TurboTrace®
Parallel
Sequential



LC/MS

Can this Handle Dirty Samples?

Typical Cartridge can have problems!

- **6ml cartridge**
 - Doesn't do well
 - Frit Surface Area is too small

Yes, A Cartridge will work

- **25ml cartridge**
 - Does well
 - 3X the Frit Surface Area



Glass/Plastic Wool

- Irregular random stranding
- Slows Particles to the Uniform Frit
- Prevents Clogging




Prepping the 6ml Cartridge with Plastic Filtration Wool

6ml cartridge with Plastic wool

- Take a little and push it into the barrel of the syringe until it touches the cartridge Frit
- The Sample will not clog, it will take longer to process




Dirty Sample from a Customer




Industrial 433 Matrix 250ml





6ml and 25ml Cartridges



250 ml run to completion on 6 ml cartridge with Plastic Wool



250ml run to completion 25ml cartridge



Automated Extraction and Concentration for PFAS/PFOS in Drinking Water



EconoTrace® PFC

Drinking Water, Serum and samples minimal particulate
Positive Pressure Pumping



LC/MS

EconoTrace for Extraction and Concentration of PFAS/PFOS in Drinking Water

Fully Automated

**Modular and expandable
from 1 to 4 Modules**

**High Throughput Runs
2 to 8 Sample Extractions in Parallel**

**Uses Positive Pressure Pumping only
for Precise delivery of
Conditioning Solvent,
Sample, Bottle Rinse and Elution Solvent**

Automatic Bottle Rinse

Delivers extract directly to SuperVap Concentrator for final blowdown



TurboTrace PFC Automated Extraction and Concentration for PFAS/PFOS in WasteWater



TurboTrace® PFC

Wastewater, Particulate Laden Samples
Drinking water Samples
Vacuum Pump for loading samples
Positive Pressure Pumping for conditioning, rinsing and elution



LC/MS

TurboTrace for Extraction and Concentration of PFAS/PFOS in Waste Water

Fully Automated

**Modular and expandable
from 1 to 6 Modules**

**High Throughput Runs Sample
Extraction in Parallel**

**Uses Positive Pressure Pumping
for Precise delivery of
Elution and Wash Solvent**

**Uses Vacuum or
Positive Pressure Pumping
to Load Samples**

Direct delivery to Concentrator for Evaporation



Automated Concentration for PFAs

- SuperVap PFC
 - 24 positions
 - 15ml Conical vials
 - Dry Bath



Concentration Functionality

- **Self Installable**
 - Video unpacking, installation and training video
- **Preprogrammed with most common temperature settings**
- **Direct to 15ml Centrifuge Tubes**
- **Dry bath heating element**
- **Time based endpoint**
- **Savable temperature log**
- **No Teflon Components**



TurboTrace Parallel Sequential PFC Automated Extraction and Concentration for PFAS/PFOS in Drinking Water and WasteWater



TurboTrace® PFC Parallel Sequential


Wastewater, Particulate Laden Samples
Drinking water Samples
Vacuum Pump or Positive Pressure for loading samples
Positive Pressure Pumping for conditioning, rinsing and elution
Direct to SuperVap Concentrator



LC/MS

EPA 1633

Turbo Trace Parallel/Sequential



Various Blanks EPA 1633

	ch.1	ch.2	ch.3	ch.4	ch.5	stopcock	sample bottle
Name	Calc. Conc.	Calc. Conc.	Calc. Conc.	Calc. Conc.	Calc. Conc.	Calc. Conc.	Calc. Conc.
ADONA	0.00	0.00	0.00	0.00	0.00	0.00	0.00
PFPeS	0.00	0.00	0.00	0.00	0.00	0.00	0.00
PFOA	0.03	0.01	0.02	0.02	0.00	0.28	0.01
PFNA	0.02	0.03	0.01	0.03	0.05	0.00	0.01
8-2 FTS	0.01	0.00	0.00	0.02	0.01	0.00	0.00
PFDA	0.02	0.02	0.03	0.01	0.02	0.04	0.00
PFDoS	0.00	0.00	0.00	0.00	0.00	0.00	0.00
PFOS	0.03	0.02	0.04	0.02	0.01	0.03	0.05
6-2 FTS	0.01	0.00	0.00	0.00	0.00	0.00	0.02
PFHxS	0.00	0.00	0.01	0.03	0.03	0.05	0.06
PFTTrDA	0.01	0.00	0.00	0.04	0.04	0.00	0.03
N-EtFOSAA	0.00	0.00	0.00	0.00	0.00	0.00	0.00
PFHpS	0.00	0.00	0.02	0.03	0.03	0.00	0.00

TurboTrace for Extraction and Concentration of PFAS/PFOS in Drinking/WasteWater

Fully Automated

**Modular and Scalable expandable
from 1 to 6 Modules**

**Run 1 to 6 samples simultaneously,
up to 30 sequentially**

**Each Module can Run 5 Samples in a
Sequence**


**Uses Positive Pressure Pumping
for Precise delivery of
Elution and Wash Solvent**

Uses Vacuum with liquid sensor to Load Dirty Samples


Direct Extract delivery to Concentrator for Evaporation



EPA 533



Semi-Automated Solid Phase Extraction front end LC/MS



EZPFC



LC/MS


Sample Analysis Work Flow

Automated Sample Prep Time

= 80 Minutes



Solid Phase Extraction
35 Minutes



Solid Phase Extraction
35 Minutes



Concentration
45 Minutes

Semi Automated Sample Prep Time

= 80 Minutes



Reasons for Semi-Automated SPE

- **Closed System**
- **12 samples in parallel**
- **Reduced Actions / Easy to Use**
- **Simplified procedures**
- **Semi-Automated versus Manual protocols = Reproducibility**
- **Increased Sample Throughput**
- **Low cost compared to Automated solutions**



Objective for Semi Automation

- **Incorporate automated system capabilities to the semi-automated platform.**
- **Our goal is to establish several Standard Operating Procedures (SOPs) for the testing lab through the utilization of a solitary extraction platform.**
- **In order to optimize time and minimize mistakes, it is crucial to decrease the amount of manual steps involved.**



Goal

- **Low Cost**
 - Extract and Concentrate a batch of 24 samples 2 to 3 hours
- **Self Installable**
 - Unpacking and Installation/training video
- **Easy to Operate**
 - No Computers or Electronics to fail or maintain
- **Semi - Automated**
 - Hyphenates the entire Solid Phase Extraction Process – Extraction, Bottle Rinse, Inline Drying and Direct to 15ml Centrifuge tube for easy Concentration
- **Fast**
 - The fastest sample processing available for SPE
 - Run up to 12 samples simultaneously
 - Vacuum for fast loading of large volume samples
 - Unattended Sample loading walkaway time
- **Closed system**
 - Eliminate potential outside contamination

Goal


- **Efficient**

- Uses all SPE cartridge sizes
- Dedicated manifold for cartridge conditioning and sample loading
- Dedicated manifold for extraction and extracts
- Separates Organic from Aqueous waste
- Vacuum cartridge drying, Nitrogen cartridge drying or combined
- Automated Bottle Rinse and Elution
- Inline Extract Drying
- Small number of components to clean

- **Low Capital Expense**



Drinking Water,
Wastewater,
Particulate laden samples



EZPFC 12 sample


System Components

No Teflon


**Tubing - High Density
Polyethylene**

Fittings – Delrin


**Cartridge Adapters –
Medical Grade
Polypropylene**



System Layout



Control Valve Layout



Automated Concentration for PFAs

- SuperVap PFC
 - 24 positions
 - 15ml Conical vials
 - Timed Endpoint




Clean up is easy with no cross contamination

- Back Flush the sample line into the original sample bottle with an IPA non-Teflon squirt bottle.
- Wash the inside of the bottle cap with IPA squirt bottle
- Wash Cartridge Adapters with IPA squirt bottle or sonicate in a beaker
- Ready for the next 12 samples

EPA 1633

EZPFC



- **EZPFC and SuperVap systems are easy to use and install**
 - Complete Water Sample Prep Workflow
- **Low cost, High throughput, Low maintenance solution**
- **EZPFC Extractions and Concentration**
 - Closed System Reduces Contamination
 - Reduces Human error



Summary

- **FMS semi-automated SPE and SuperVap systems deliver consistent, reproducible results**
- **Handles a wide range of Sample sizes and matrix types**
- **Uses all SPE Cartridge sizes**
- **Comply with existing methods that require vacuum, positive pressure and precise delivery of sample and solvents**



Automated SPE for PFAS

- **EconoTrace**
 - Fully automated, positive pressure system, automated bottle rinse, direct delivery to SuperVap Concentrator
 - Designed for cleaner samples
- **TurboTrace**
 - Fully automated, Vacuum pump for loading “dirty samples”, automated bottle rinse,
 - Direct delivery to SuperVap Concentrator

Semi Automated SPE for PFAS

- **EZPFC**
 - 2 EZPFC 12 position systems and the SuperVap PFC Concentrator for 24 samples
 - Run a batch of 24 in 2 hours

