

Low Cost, Compact, High Throughput, Automated Sample Concentration and Evaporation

SuperVap®




Agenda

- Company Overview
- SuperVap Overview
- Questions



- Fluid Management Systems
 - Founded in 1986
 - Focus
 - Automating the Sample Prep Process
 - Markets
 - Agricultural
 - Environmental
 - Clinical
 - Food and Beverage
 - Pharmaceutical
 - Petrochemical

Made in the USA



Solid Phase Extraction



Pressurized Liquid Extraction



Automated Sample Cleanup



Direct to Vial Concentration



Sample Handling

- Large Volume Concentration
- Concentrate/Evaporate up to 6 Samples
- Sample Sizes up to 220ml
- Automatic Endpoint Detection and Nitrogen Shutoff for each Vessel
- Compact Size



Concentration Vessels



Direct to GC vial Vessel



Sample Handling

- Small Volume Concentration
- Concentrate/Evaporate up to 12 Samples
- Sample Sizes up to 50ml
- Automatic Endpoint Detection and Nitrogen Shutoff for each Vessel
- Compact Size



Concentration Vessels



Sample Handling

- Small Volume Evaporation
- Evaporate up to 24 Samples
 - Sample Size Format
 - 2ml vial
 - 4ml vial
- Evaporate up to 12 Samples
 - Sample Size Format
 - 20ml Vial
 - ASE 40ml Vial
 - ASE 60ml Vial
- Timed or Manual Nitrogen Shutoff for each Vessel
- Compact Size



Evaporation Vessels



Automated Concentration for PFAs



Automated Concentration for PFAs

- SuperVap PFC
 - 24 positions
 - 15ml Conical vials



Automated Concentration for PFAs

- SuperVap PFC
 - 12 positions
 - 50 ml Conical vials



No Waterbath


- Dry Heating Assembly
 - Robust Endpoint Sensors
 - No Water dripping into the vessel as in a water bath




Easy to Use

- Touch Screen Programming
 - Programmable Temperature
 - 0° to 100° C depending upon model
 - End point liquid level sensor
 - Sensor sees the liquid has reached desired volume turns off Nitrogen
 - Timed End point
 - Set the time for the Nitrogen and Heat to turn off


Front View



Rear View



New Home Screen



File Storage

Open / Create Evap File

7 Files on disk. Max=200

T_55.EVP	
T_60.EVP	
T_65.EVP	
T_70.EVP	
TEST_N2_6.EVP	
T_45_DRY.EVP	
T_50.EVP	

A	B	C	D	E	F	G	H
I	J	K	L	M	N	O	P
Q	R	S	T	U	V	W	X
Y	Z	_	0	1	2	3	4
5	6	7	8	9		BK	
Create						CLR	

Delete Select Cancel

Selecting a File

Open / Create Evap File


7 Files on disk. Max=200

T_55.EVP	
T_60.EVP	
T_65.EVP	
T_70.EVP	
TEST_N2_6.EVP	
T_45_DRY.EVP	
T_50.EVP	


A	B	C	D	E	F	G	H
I	J	K	L	M	N	O	P
Q	R	S	T	U	V	W	X
Y	Z	_	0	1	2	3	4
5	6	7	8	9		BK	
Create						CLR	

Delete Select Cancel

Running a Method



Running a Method



Vessels to Completion

File **T_50.EVP** **Open** **X**

Stp	Time	Tmax	Desc
01	5.2	46	PREHEAT
02	.1	50	EVAP

Stp **2** Temp **52.04** Elapsed **16.600**

Mode
☐ Time
☒ Sensor

☒ Heat

☐ No Alarm

Tubes
Mode: ☒ **AUTO** ☐ **MANUAL**

☐ 1 ☐ 2 ☐ 3
☐ 4 ☐ 5 ☐ 6

Run **Stop**

Halt **Offset** **Single** **STOP**


Plot **Exit**

Manual Switch
☐ 1 ☒ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6
OFF **1,3,**

Temperature Plot



Temperature at any point



Manual Operation



Create a DataLog to Save a Plot

Remaining Time: DataLog

Hours Minutes Seconds

☐ File (Sample) Name:


☐ Log Frequency (x4)

Sec

A	B	C	D	E	F	G	H
I	J	K	L	M	N	O	P
Q	R	S	T	U	V	W	X
Y	Z	_	0	1	2	3	4
5	6	7	8	9	<input type="button" value="BK"/>		
					<input type="button" value="CLR"/>		

Status:

Create or Open a File



Programming a Method

The screenshot displays the FMS programming interface. At the top, a 'File' menu shows 'T_50.EVP' with 'Open' and 'X' buttons. Below this is a table with columns 'Stp', 'Time', 'Tmax', and 'Desc'. The table contains three rows: '01 5.2 46 PREHEAT', '02 .1 50 EVAP', and '03 .1 50 EVAP'. The third row is highlighted. To the right of the table is a 'Tubes' section with six numbered icons (1-6) arranged in two rows of three. At the bottom, there are buttons for 'Copy', 'Delete', 'Edit', 'Save', and 'Exit'. A 'Mode' section at the bottom right shows 'Time' and 'Sensor' options with corresponding icons, and a 'Heat' option with a red square icon.

Stp	Time	Tmax	Desc
01	5.2	46	PREHEAT
02	.1	50	EVAP
03	.1	50	EVAP

Tubes

1 2 3
4 5 6

Mode

Time
Sensor
Heat

Saving a Method

Save Evap File

Save

Save As...


Cancel

System Configuration

Configurations

Reboot	EVP6 Unit Name/ID
Clock	LCD Timeout 11
<input type="checkbox"/> Use as Slave	Exit

Vial Concentrator Interface



Vial/Centrifuge Interface

File

Stp	Time	Tmax	Desc
01	20	60	PREHEAT
02	9999	60	EVAP

Stp Temp Elapsed

☐ Heat

Vials
Mode: ☒ AUTO ☐ MANUAL

☐ 1 ☐ 2 ☐ 3 ☐ 4

System Calibration


Temperature Calibration

Select Cal Point ☐ 1 ☐ 2

Enter measured temperatures deg

deg
 mv

Automated Extraction and Concentration




Alkane Recovery Data

Compound	Percent Recovery
Nonane (C9)	75%
Decane (C10)	77%
Dodecane (C12)	88%
Tetradecane (C14)	92%
Hexadecane (C16)	95%
Octadecane (C18)	97%
Nonadecane (C19)	97%
Eicosane (C20)	98%
Docosane (C22)	98%
Tetracosane (C24)	99%
Hexacosane (C26)	98%
Octacosane (C28)	97%
Triacontane (C30)	96%
Hexatriacontane (C36)	97%


Typical PAH Recovery Data

Compound	Percent Recovery
Naphthalene	78%
2-Methylnaphthalene	102%
Acenaphthylene	83%
Acenaphthene	83%
Fluorene	87%
Phenanthrene	89%
Anthracene	89%
Fluoranthene	93%
Pyrene	90%
Benzo[a]anthracene	86%
Chrysene	95%
Benzo[b]fluoranthene	90%
Benzo[k]fluoranthene	93%
Benzo[a]pyrene	89%
Indeno[1,2,3-cd]pyrene	90%
Dibenzo[a,h]anthracene	89%
Benzo[g,h,i]perylene	91%


EPA 537.1



PFAS Background



EPA 533



EPA 1633 SuperVap 12 PFC

Name	Average (%)	RSD (%)	EPA PASSING
			Window
13C2-4-2 FTSA	94.7	6.5	40-200
13C2-6-2 FTS	95.6	7.3	40-200
13C2-8-2 FTSA	85.0	7.0	40-300
13C2-PFDoDA	104.0	7.0	10-130
13C2-PFTDA	100.1	3.5	10-130
13C3-HFPO-DA	101.4	4.0	40-130
13C3-PFBS	84.7	5.9	40-135
13C3-PFHxS	93.5	7.1	40-130
13C4-PFBA	90.5	7.7	5-130
13C4-PFHpA	89.6	8.7	40-130
13C5-PFHxA	94.4	7.4	40-130
13C5-PFPeA	108.2	3.4	40-130
13C6-PFDA	82.4	8.3	40-130
13C7-PFUnA	92.2	6.8	30-130
13C8-PFOA	99.4	8.3	50-200
13C8-PFOS	93.3	7.1	50-200
13C8-PFOSA	97.3	7.2	40-130
13C9-PFNA	83.9	9.2	40-130
2H3-N-MeFOSA	99.2	3.7	10-130
2H3-N-MeFOSAA	99.1	4.9	40-170
2H5-N-EtFOSA	97.6	9.1	10-130
2H5-N-EtFOSAA	91.4	6.3	25-135
2H7-MeFOSE	95.1	4.1	10-130
2H9-EtFOSE	97.2	6.3	10-130

The SuperVap in Summary

- **Waterless Bath**
 - Dry Heater Assembly
 - Reduce Errors with Robust Sensors
 - No Water dripping into the vessel from condensation
 - Eliminates bacteria, biofilms, and cross-contamination from water
- **Fully Automated**
- **Compatible with all FMS automated systems**
 - Standalone or Integrated
 - Fully automated Sample Prep WorkFlow



The SuperVap in Summary

- Sample Handling
 - Direct to Vial
 - Reduces error associated with transfer steps
 - Concentrate/Evaporate up to 6 Samples
 - Sample Sizes up to 220ml
 - Concentrate/Evaporate up to 12 Samples
 - Sample sizes up to 50ml
 - Automatic Endpoint Detection and Nitrogen Shutoff for each Vessel
 - Timed Endpoint



The SuperVap in Summary

- Variety of Models and Vessels
 - Direct to GC/LC vial
 - 500ul
 - 1ml low volume GC/LC vial
 - 2ml GC/LC vial
 - 4ml GC/LC vial
 - 20ml VOA vial
 - ASE 40ml and 60ml VOA vials
 - 15ml Centrifuge Tube
 - 50ml Centrifuge Tube

The SuperVap in Summary

- Easy to Use Touch Screen
- Smallest Footprint
- Inexpensive Glassware
- Consistent Reproducible Results
- Reduces Errors



Questions

